


Product datasheet for PS100063

pCMV6-Entry2 Mammalian Expression Vector

Product data:

Product Type: Mammalian Expression Vectors
Description: pCMV6-Entry2, mammalian vector with C-terminal Myc-DDK tag, 10ug

Product images:



pCMV6-Entry2

										Kozac Consensus					
<i>EcoR I</i>		<i>BamH I</i>		<i>Kpn I</i>		RBS		<i>Sgf I</i>		<i>Asc I</i>		<i>Bgl II</i>			
CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGCGCCAGATCT															
<i>Hind III</i>		<i>Nhe I</i>		<i>Rsr II</i>		<i>Mlu I</i>		<i>Not I</i>		<i>Xho I</i>		Myc.Tag			
CAAGCTTAAGTACTAGCTAGCGGACCG <u>ACG CGT</u> <u>ACG CGG CCG CTC</u> <u>GAG CAG AAA</u> <u>CTC ATC TCA GAA GAG</u>															
T R T R P L E Q K L I S E E															
				<i>EcoR V</i>				DDK.Tag				<i>Pme I</i>		<i>Fse I</i>	
GAT CTG GCA GCA AAT <u>GAT ATC</u> <u>CTG GAT TAC AAG GAT GAC GAC GAT AAG</u> <u>GTT TAA</u> <u>ACGGCCGGCC</u>															
D L A A N D I L D Y K D D D D K V Stop															



[View online >](#)